

Headrest / Neckrest for HD Balance



There are two types of support for the patient's head / neck, Headrest (A) and Neckrest (B).





Headrest

The headrest gives a slightly larger support for the head and the shape provides a comfortable support to the sides. The lower part is slightly bevelled in order to not give pressure in the nape. The headrest can also be turned upside down if a thicker and rounder shape in the nape is desirable.

The headrest is available in three models:

Basic: With a headrest cushion without side adjustability (C)

Side: With a side adjustable headrest cushion (D)

Side+: With side adjustable headrest cushion and side support brackets (right and left) (E)



Neckrest

The neckrest (B) is slightly smaller and acts as a support in the nape of the more active user.

The neckrest is available in two models:

Side: With a side adjustable neckrest cushion (F)

Side+: With a side adjustable neckrest cushion and side support brackets (right and left) (G)



Side Support

Side support is used when the patient has an increased need for lateral stability on one or both sides. The side supports are reversible and can be used on the left or right side. Side support cushions are the same for the headrest and neckrest. The bends, however, differ in the form.


The side support is available in two models: **Fix:** Side support with fixed angle for headrest (H) **Vario:** Side support with adjustable angle for headrest (I) and for necksupport (J)






Adjustments

The head- / neckrest consists of two parts; an extension bar (K) which is monted in the bracket on the back of the wheelchair , and a head- / neckrest bar (L) with a cushion. The head- / neckrest can be adjusted in height (24 cm) – depth (10 cm) using knobs. The angle of the cushion (M) and the angle between the headrest bar and extension bar are adjustable without tools.







!qiT

Use the profile stop that enables to always set the head- / neckrest in the same position (accessory).



Dimensions - cushions:

Table1 - Dimensions - cushions			
Model	height [cm]	width [cm]	thickness [cm]
Headrest	16	31	5 (centre of the back)
Neckrest	13	28	6 (centre of the back)
Side support	13	11,5	2,5

Keep in mind:

It is possible to provide the head-/ neckrest with a longer head-/ neckrest bar (O) to increase the adjustability forward (10 cm).

This is a special design which needs to be evaluated in terms of risks. Evaluate if the benefits of this special design overcome the additional risks.



Covers

Head- and neckrest comes as standard with grey Trevira CS fabric upholstery. They can also be ordered in black Trevira CS or black disinfactable Dartex. All covers are removable and washable at 60° C. Covers can also be ordered separately.

Article numbers

All article numbers ars found in the Prescription Form, art nr 95750-1.