

Wedge side supports HD Balance

Art. no.: 63440

User manual / Assembly instructions

Page 1 (1)

Parts included

Wedge side supports Wide



Wedge side supports Narrow



Assembly



1. Remove the current side supports by loosening the three screws and pulling the side support off the armrest.
2. Mount the wedge side supports as seen in the picture. Tighten the screws.
3. The picture shows the wedge side supports mounted on the chair.

Care instructions

The wedge side supports can be cleaned by wiping with soap and water. For further information see Maintenance Instructions, art. no. 95730-1.