

Date		Order number	
Customer		Contact person	
Customer tel.		Contact person tel.	
Customer email		Contact person email	
User		Contact at HD Rehab	

- This Prescription Form is used when ordering a **single** wheelchair (with accessories) in a particular configuration.
- All parts that are included in standard configuration are marked with a bold lined tickbox (☐). These will be delivered with the wheelchair if no other option is selected.
- For further information on different configurations and accessories, see the Prescriber Information (doc.nr. 95715-1) and each specific Info Sheet. User manuals and assembly instructions are also available at hdrehab.com.

1. CHASSIS

Info Sheet 95937-1 Chassis

	Seat width		Standard		Reinforced
HD Balance 24	38	☐	60210-6	☐	60220-6
Self propelled version	42	☐	60310-6	☐	60320-6
Seat depth 43-52cm	46	☐	60410-6	☐	60420-6
(with seat cushion)	50	☐	60510-6	☐	60520-6
HD Balance 24	38	☐	60215-6	☐	60225-6
Self propelled version	42	☐	60315-6	☐	60325-6
Extended seat depth 48-57cm	46	☐	60415-6	☐	60425-6
(with seat cushion)	50	☐	60515-6	☐	60525-6
HD Balance 16	38	☐	60250-6	☐	60260-6
Transport model	42	☐	60350-6	☐	60360-6
Seat depth 43-52cm	46	☐	60450-6	☐	60460-6
(with seat cushion)	50	☐	60550-6	☐	60560-6
HD Balance 16	38	☐	60255-6	☐	60265-6
Transport model	42	☐	60355-6	☐	60365-6
Extended seat depth 48-57cm	46	☐	60455-6	☐	60465-6
(with seat cushion)	50	☐	60555-6	☐	60565-6



2. WHEELS & CASTORS

Info Sheet 95938-1 *Wheels*

		HD Balance 24	HD Balance 16
PUR (Polyurethane) drive wheels <i>Solid tyres with aluminium hand rims</i>	☒	67505-1	67551-1
Flexel drive wheels <i>Solid tyres with aluminium hand rims</i>	☐	67506-1	67552-1
Pneumatic drive wheels <i>With aluminium hand rims</i>	☐	67507-1	67553-1




Castors 175 x 43mm	☒	67310
Castors 175 x 28mm	☐	67316
Castors 125 x 27mm, extended fork	☐	67318
Castors Flexel 200 x 50mm <i>Conflict with footrest in 90° position</i>	☐	67314
Castors 200 x 28mm <i>Conflict with footrest in 90° position</i>	☐	67317



3. HAND RIMS (HD Balance 24)

Info Sheet 95938-1 *Wheels*

Aluminium	☒	15491
Stainless steel	☐	15492
Soft grip, plastic coated	☐	15493



4. PUSH BAR / HANDLES

Info Sheet 95939-1 *Push bar / Handles / Caregiver brakes*

Side unit, standard, left and right <i>With caregiver brakes</i>	☒	62330-1, 62330-2
Side unit, extended +12cm, left and right <i>With caregiver brakes. Longer wires for brakes, seat tilt, and back recline included.</i>	☐	62330-5, 62330-6
Side unit, angle adjustable, left and right <i>With caregiver brakes</i>	☐	62392-1, 62392-2



		sw38	sw42	sw46	sw50
Push bar, standard	☒	62341-1	63341-1	64341-1	65341-1
Push bar, extra depth (+8cm)	☐	62342-1	63342-1	64342-1	65342-1
Push bar, extended & bent <i>Not in combination with Side units extended or angle adjustable</i>	☐	62343-1	63343-1	64343-1	65343-1
Separate operating handles <i>Not in combination with Side units angle adjustable</i>	☐	10295			



5. SEAT CUSHION

Info Sheet 95940-1 *Seat cushions*

		sw38	sw42	sw46	sw50
Seat cushion standard, Grey Trevira	☐	68201	68301	68401	68501
Seat cushion standard, Black Trevira	☐	68203	68303	68403	68503
Seat cushion standard, Hygiene	☐	68204	68304	68404	68504



Seat cushion long, +5cm, Grey Trevira	☐	68205	68305	68405	68505
Seat cushion long, +5cm, Black Trevira	☐	68207	68307	68407	68507
Seat cushion long, +5cm, Hygiene	☐	68208	68308	68408	68508

(for chassis with extended seat depth)

Elastic seat cover, Incontinence	☐	pcs	49206	49306	49406	49506
Elastic seat cover +5cm, Incontinence	☐	pcs	49208	49308	49408	49508

Elastic seat covers are fitted over existing cover. All elastic seat covers are Black Trevira.


Without seat cushion	☐	68200-8
----------------------	--------------------------	---------



6. BACKREST SYSTEM

Info Sheet 95941-1 *Backrest systems*

		sw38	sw42	sw46	sw50
Flexi-back standard	☐	62250	63250	64250	65250
Flexi-back long (+10cm)	☐	62255	63255	64255	65255
Solid back standard	☐	62260	63260	64260	65260
Solid back long (+10cm)	☐	62265	63265	64265	65265



Note that increased or decreased width for the backrest system can be selected in section 15.

7. BACK CUSHION

Info Sheet 95941-1 *Backrest systems*

		sw38	sw42	sw46	sw50	
Back cushion standard, Grey Trevira	☒	68281	68381	68481	68581	
Back cushion standard, Black Trevira	☐	68291	68391	68491	68591	
Back cushion standard, Dartex <i>with side supporting wings</i>	☐	68282	68382	68482	68582	
Back cushion flat, Grey Trevira	☐	68285	68385	68485	68585	
Back cushion flat, Black Trevira	☐	68295	68395	68495	68595	
Back cushion flat, Dartex <i>without side supporting wings</i>	☐	68286	68386	68486	68586	
Back cushion std long (+10cm), Grey	☐	68283	68383	68483	68583	
Back cushion std long (+10cm), Black	☐	68293	68393	68493	68593	
Back cushion std long (+10cm), Dartex <i>with side supporting wings</i>	☐	68284	68384	68484	68584	
Back cushion flat long (+10cm), Grey	☐	68287	68387	68487	68587	
Back cushion flat long (+10cm), Black	☐	68297	68397	68497	68597	
Back cushion flat long (+10cm), Dartex <i>without side supporting wings</i>	☐	68288	68388	68488	68588	
Without back cushion	☐	68250-8				

8. LEG RESTS

Info Sheet 95933-1 Leg rests


	Left	Right
Angle adjustable leg rest, standard	☐ 62725	☐ 62720
Angle adjustable leg rest, short	☐ 62721	☐ 62722
Angle adjustable leg rest, long	☐ 62723	☐ 62724
Disengage mechanism - knee joint <i>One sided adjustment of leg rest angle for one piece footrest / calf support / foot box</i>	☐ 62733-5	<i>RH version on request</i>
Amputee leg rest <i>Angle adjustable, including calf support</i>	☐ sw38 - 62771 sw42 - 63771 sw46 - 64771 sw50 - 65771	☐ sw38 - 62772 sw42 - 63772 sw46 - 64772 sw50 - 65772
Fixed angle leg rest, 90°, standard	☐ 63758-1	☐ 63758-2
Fixed angle leg rest, 90°, short	☐ 63758-5	☐ 63758-6
Fixed angle leg rest, 90°, long	☐ 63758-3	☐ 63758-4
Fixed angle leg rest, 105°, standard	☐ 63759-1	☐ 63759-2
Fixed angle leg rest, 105°, short	☐ 63759-5	☐ 63759-6
Fixed angle leg rest, 105°, long	☐ 63759-3	☐ 63759-4
Without leg rest <i>Calf supports and Footrests then also excluded</i>	☐ 62701-8	☐ 62702-8



9. CALF SUPPORTS



Info Sheet 95934-1 Calf supports

	sw38	sw42	sw46	sw50
Calf support separate, left	☐ 62801	☐ 63801	☐ 64801	☐ 65801
Calf support separate, right	☐ 62802	☐ 63802	☐ 64802	☐ 65802
Calf support sep., +30mm extra depth, left	☐ 62803	☐ 63803	☐ 64803	☐ 65803
Calf support sep., +30mm extra depth, right	☐ 62804	☐ 63804	☐ 64804	☐ 65804
Calf support sep., +50mm extra depth, left	☐ 62805	☐ 63805	☐ 64805	☐ 65805
Calf support sep., +50mm extra depth, right	☐ 62806	☐ 63806	☐ 64806	☐ 65806
One piece calf support	☐ 62811	☐ 63811	☐ 64811	☐ 65811
One piece calf support, +30mm extra depth <i>Order disengage mechanism for knee joint separately</i>	☐ 62812	☐ 63812	☐ 64812	☐ 65812
Calf support panel	☐ — 63891 —	— 65891 —		
Calf support strap	☐ — 63892 —			
Mounts for calf panel / calf strap, 2 pcs	☐ — 63895 —			
Without calf support left	☐ 63801-8			
Without calf support right	☐ 63802-8			










10. FOOTRESTS

Info Sheet 95935-1 Footrests

		sw38	sw42	sw46	sw50	
Separate footrest, standard, left	☐	62909	63909	64909	65909	
Separate footrest, standard, right	☐	62908	63908	64908	65908	
Separate footrest, short, left	☐	62903	63903	64903	65903	
Separate footrest, short, right	☐	62904	63904	64904	65904	
<i>For use with short leg rests</i>						
One-piece footrest, standard	☐	62951	63951	64951	65951	
<i>Order disengage mechanism for knee joint separately</i>						
One-piece footrest, short	☐	62953	63953	64953	65953	
<i>For use with short leg rests</i>						
<i>Order disengage mechanism for knee joint separately</i>						
Footbox, standard	☐	62981	63981	64981	65981	
<i>Order disengage mechanism for knee joint separately</i>						
Footbox, short	☐	62982	63982	64982	65982	
<i>For use with short leg rests</i>						
<i>Order disengage mechanism for knee joint separately</i>						
Padding, separate	☐	62991-1		64991-1		
<i>Kit with 2 pieces, for separate footrests</i>						
Padding, one-piece	☐	62993-1	63993-1	64993-1	65993-1	
<i>For one-piece footrest or footbox</i>						
Footrest side support, left	☐	65987				
Footrest side support, right	☐	65986				
<i>For use with standard length footrests</i>						
Mounting of footrest side support	☐	65987-8				
Heel strap, kit with 2 pieces	☐	62995-1				
<i>For use with separate footrests</i>						
Mounting of heel straps	☐	62995-8				
Without footrest left	☐	63909-8				
Without footrest right	☐	63908-8				

11. ARMRESTS

Info Sheet 95936-1 Armrests

	Left	Right		
Armrest standard	☒ 62401	☒ 62402		
Armrest wide	☐ 63401	☐ 63402		
Armrest standard low	☐ 62405	☐ 62406		
Armrest wide low	☐ 63405	☐ 63406		
Armrest standard with wedge side support 7°	☐ 62411	☐ 62412		
Armrest wide with wedge side support 7°	☐ 63411	☐ 63412		
Armrest standard with wedge side support 4°	☐ 62413	☐ 62414		
Armrest wide with wedge side support 4°	☐ 63413	☐ 63414		
Armrest standard with long side support	☐ 62415	☐ 62416		
Armrest wide with long side support	☐ 63415	☐ 63416		
Hemiplegia armrest, high	☐ 63491	☐ 63492		
Hemiplegia armrest, low	☐ 63493	☐ 63494		
Side support standard with post, w/o arm pad	☐ 62451-1	☐ 62451-2		
Post with cap, without arm pad	☐ 62453-1	☐ 62453-1		
Wedge side support 7°	☐ 63463	☐ 63464		
Wedge side support 4°	☐ 63461	☐ 63462		
Long side support	☐ 63465	☐ 63466		
Coaming for armrest, standard	☐ 62471	☐ 62472		
Coaming for armrest, wide	☐ 63471	☐ 63472		
Mounting of armrest coaming	☐ 62471-8			
Armrest locks, pair <i>Locks the armrests to the seat frame</i>	☐ 62459			
Mounting of armrest locks	☐ 62459-8			
Without armrest	☐ 62400-8	☐ 62400-8		

ALTERNATE CONFIGURATIONS

Alternate configurations of HD Balance wheelchairs must be evaluated for risk. Evaluate if there is an increased risk of tipping and/or pinching for the particular user with the desired modification. If there is increased risk, make an informed decision as to whether the alternate configuration is still warranted. **Inform everyone involved about the risks.**

12. BALANCE POINT

Info Sheet 95942-1 Balance Point

Forward mounted seat unit	☐	67056-8	See 95758-1 Balancing modes, 95759-1 Reference table
Back mounted seat unit	☐	67057-8	balancing configurations, and 95707-1 Risk information
Drive wheels rear position (5 cm back)	☐	67551-8	

13. SEAT HEIGHT AND SEAT DEPTH

Info Sheet 95943-1 Seat height, 95937-1 Chassis

Increased seat height +5 cm	☐	67052-8	See 95758-1 Balancing modes, 95759-1 Reference table
Lowered seat height -3 cm	☐	67053-8	balancing configurations, and 95707-1 Risk information
20" wheels for HD Balance 24 <i>Lowers seat height approx. 4 cm</i>	PUR ☐	67515-2	
	Flexel ☐	67516-2	
	Pneumatic ☐	67518-2	
125x27 castors for HD Balance 24 <i>Only used in combination with 20" wheels</i>	☐	67315	
Back spacers, kit with 4 pieces <i>Can be mounted on one side if desired</i>	☐	69791-5	See 95707-1 Risk information 
Mounting of back spacers	☐	69791-8	

14. SEAT TILT AND BACK RECLINE

Info Sheet 95944-1 Seat tilt and Back recline

See 95707-1 Risk information. Standard configuration has 0° - 20° seat tilt and 90° - 120° back recline.

Seat tilt -5° - 20° <i>Allows 5° forward tilt</i>	☐	67820-8	
Increased seat tilt 0° - 30° <i>Note that drive wheels should be in rear position</i>	☐	67823-8	
Increased seat tilt -5° - 30° <i>Allows 5° forward tilt</i> <i>Note that drive wheels should be in rear position</i>	☐	67824-8	
Increased back recline 90° - 135° <i>Note that drive wheels should be in rear position</i>	☐	67705-8	
Lockable back recline angle	☐	62369	
Fixed back recline angle <i>Control and gas spring to reset are not included</i>	☐	Recline angle on delivery (90° - 105° possible):	☐

15. INCREASED OR DECREASED BACK WIDTH

Info Sheet 95945-1 *Back System - increased or decreased width*

		sw38	sw42	sw46	sw50
Wider back system +4 cm, incl. wider back cushion	☐	63298-8	64298-8	65298-8	On request
Wider back system +8 cm, incl. wider back cushion	☐	64299-8	65299-8	On request	On request
<i>Causes certain limitations for the backrest height setting and armrest depth</i>					
Narrower back system -4 cm, incl. narrower cushion	☐	On request	62297-8	63297-8	64297-8

16. CONFIGURATION FOR SEVERE LOAD

Info Sheet 95946-1 *Configurations for severe load from user*

Shock absorbing back ☐

Available for both adjustable and fixed back recline angle

17. BRAKES, SPECIAL VERSIONS

Info Sheet 95947-1 *Brakes, Special versions*

User brakes for HD Balance 16 ☐ 67130-1

User brakes with extended lever ☐ 67115-5

For HD Balance 24

18. OUTER DIMENSIONS

Info Sheet 95938-1 *Wheels, 95719-1 Measures*

Castor pin narrow ☐ 69991-1

Reduces total width (in combination with narrow position for hand rims)

Narrow position for hand rims ☐ 15499-8

Reduces total width, not possible for user to drive

ACCESSORIES

19. HEADREST / NECK REST


Info Sheet 95932-1 Headrest / Neck rest

	Grey Trevira	Black Trevira	Dartex
Headrest Basic	☐ 69108-1	☐ 69108-3	☐ 69108-2
Headrest Side, sideways adj.	☐ 69101-1	☐ 69101-3	☐ 69101-2
Headrest Side+ <i>With mounts for side supports Side supports ordered separately</i>	☐ 69104-1	☐ 69104-3	☐ 69104-2
Side support for headrest, Vario <i>Adjustable angle</i>	☐ pcs 69182-1	☐ pcs 69182-3	☐ pcs 69182-2
Side support for headrest, Fix <i>Fixed angle</i>	☐ pcs 69181-1	☐ pcs 69181-3	☐ pcs 69181-2
Headrest cover <i>Black Trevira</i>		☐ pcs 49911	
Neck rest Side	☐ 69111-1	☐ 69111-3	☐ 69111-2
Neck rest Side+ <i>With mounts for side supports Side supports ordered separately</i>	☐ 69114-1	☐ 69114-3	☐ 69114-2
Side support for neck rest, Vario <i>Adjustable angle</i>	☐ pcs 69191-1	☐ pcs 69191-3	☐ pcs 69191-2
Special variant: mounting bar extended +10cm		☐ 12780-1	



20. THORACIC SUPPORT & SHOULDER SUPPORT


Info Sheet 95931-1 Thoracic support, 95954-1 Shoulder support

		Grey Trevira	Black Trevira	Dartex	
Thoracic supp. Standard, narrow pos.		☐ pcs 69301-1	☐ pcs 69301-3	☐ pcs 69301-2	
Thoracic supp. Standard, wide pos.		☐ pcs 69302-1	☐ pcs 69302-3	☐ pcs 69302-2	
Thoracic support Adjustable Short		☐ pcs 69311-1	☐ pcs 69311-3	☐ pcs 69311-2	
Thoracic support Adjustable Long		☐ pcs 69312-1	☐ pcs 69312-3	☐ pcs 69312-2	
Thoracic support Adjustable Mini			☐ pcs 69332-3		
Thoracic support Standard, narrow position, with Swing	left right	☐ pcs 69305-1 ☐ pcs 69306-1	☐ pcs 69305-3 ☐ pcs 69306-3	☐ pcs 69305-2 ☐ pcs 69306-2	
Thoracic support Standard, wide position, with Swing	left right	☐ pcs 69307-1 ☐ pcs 69308-1	☐ pcs 69307-3 ☐ pcs 69308-3	☐ pcs 69307-2 ☐ pcs 69308-2	
Thoracic supp. Adj. Short, narrow position, with Swing	left right	☐ pcs 69315-1 ☐ pcs 69316-1	☐ pcs 69315-3 ☐ pcs 69316-3	☐ pcs 69315-2 ☐ pcs 69316-2	
Thoracic supp. Adj. Short, wide position, with Swing	left right	☐ pcs 69317-1 ☐ pcs 69318-1	☐ pcs 69317-3 ☐ pcs 69318-3	☐ pcs 69317-2 ☐ pcs 69318-2	
Thoracic supp. Adj. Long, narrow position, with Swing	left right	☐ pcs 69325-1 ☐ pcs 69326-1	☐ pcs 69325-3 ☐ pcs 69326-3	☐ pcs 69325-2 ☐ pcs 69326-2	
Thoracic supp. Adj. Long, wide position, with Swing	left right	☐ pcs 69327-1 ☐ pcs 69328-1	☐ pcs 69327-3 ☐ pcs 69328-3	☐ pcs 69327-2 ☐ pcs 69328-2	
Thoracic supp. Adj. Mini, narrow position, with Swing			left ☐ pcs 69323-1 right ☐ pcs 69323-2		
Thoracic support mount		☐ pcs 69351-1			
Mounting of thoracic support mount		☐ 69351-8			
Shoulder support, Swing, complete with mount			left ☐ 69201-3 right ☐ 69202-3		
Mounting of shoulder support		☐ 69201-8			

21. POMMELS

Info Sheet 95948-1 *Pommels*

Leg divider, small	☐	69605-2
Leg divider, large	☐	69606-2
Knee divider, small	☐	69615-2
Knee divider, large	☐	69616-2
Pommel mount	☐	69620-1
Mounting of pommel mount	☐	69620-8



22. TRAY TABLE

Info Sheet 95930-1 *Tray tables*

		sw38	sw42	sw46	sw50
Tray table, PC (smoke coloured)	☐	69521	69531	69541	69551
Tray table, PC (smoke), with pads	☐	69522	69532	69542	69552
Tray table, Wood (oak)	☐	69525	69535	69545	69555
Tray table, Wood (oak), with pads	☐	69526	69536	69546	69556
Support table, padded	☐	69596		69598	
Mug tray (oak)	left ☐	69581	right ☐	69582	
Tray table mounts, pair	☐	69515-1			
Tray table mounts w/ locking, pair	☐	69515-5			
Mug tray mount w/ locking	left ☐	69585-1	right ☐	69585-2	
Mounting of tray mounts	☐	69515-8			



23. POSITIONING BELT & MOUNTS

Info Sheet 95906-1 *Pos. belt*, 95949-1 *Belt mounts*, 95959-1 *Harness mounts*

2-point Positioning belt
Order 1 set of belt mounts separately

☐ 13545



4-point Positioning belt
Order 2 sets of belt mounts separately

☐ 13546



Belt mounts, set with 2 pieces

☐ sets 69710-1



Mounting of belt mounts

☐ 69710-8

Harness mounts for Flexi-back, pair
Including strap holders for lower straps

☐ 69720-5



Harness mounts for Solid back, pair
Including strap holders for lower straps

☐ 69730-5



Mounting of harness mounts

☐ 69720-8 (Flexi)
69730-8 (Solid)

24. SPOKE PROTECTION & WHEEL COVERS


Info Sheet 95938-1 *Wheels*

Spoke protection for HD Balance 24, pair

☐ (24" - 15471, 20" - 15472)

Mounting of spoke protection

☐ 15475-8



Wheel covers for HD Balance 16, pair

☐ 15690

Mounting of wheel covers

☐ 15690-8



25. ELECTRICAL FUNCTIONS

Info Sheet 95951-1 *Electrical functions*

Electric seat tilt and back recline
*Only in combination with standard seat tilt (0-20°)
and back recline (90-120°)*

☐ 69801




Electric seat tilt
Only in combination with standard seat tilt (0-20°)

☐ 69802

26. OTHER ACCESSORIES

Side pads 1 cm, pair <i>Reduces seat width 1 cm per side. Info Sheet 95937-1.</i>	☐	69680
Spacers for side supports 1 cm, pair <i>Reduces seat width 1 cm per side. Info Sheet 95937-1.</i>	☐	62458-1
Side support padding, pair <i>Soft surface against thigh. Info Sheet 95936-1.</i>	☐	69685
Back wedge <i>Pushes back cushion forward on one side. Info Sheet 95958-1.</i>	☐	69698-1
Profile stop <i>Assists in mounting supports (head/neck/thoracic) in correct position. Info Sheet 95932.</i>	☐	pcs 12799
Bag hook <i>Info Sheet 95957-1.</i>	☐	pcs 13605
Mounting of bag hook	☐	13605-8
IV pole <i>Info Sheet 95956-1 and Risk Information 95707-1.</i>	☐	69760
IV pole mount	left ☐	69755-1
	right ☐	69755-2
Mounting of IV pole mount	☐	69750-8
Oxygen tank holder 2 litre <i>Info Sheet 95955-1 and Risk Information 95707-1.</i>	☐	56905
Mounting of oxygen tank holder	☐	56904-8



27. NOTES

