

| | | | |
|----------------|----------------------|----------------------|----------------------|
| Date | <input type="text"/> | Order number | <input type="text"/> |
| Customer | <input type="text"/> | Contact person | <input type="text"/> |
| Customer tel. | <input type="text"/> | Contact person tel. | <input type="text"/> |
| Customer email | <input type="text"/> | Contact person email | <input type="text"/> |
| User | <input type="text"/> | Contact at HD Rehab | <input type="text"/> |

- This Prescription Form is used when ordering a **single** wheelchair (with accessories) in a particular configuration.
- All parts that are included in standard configuration are marked with a bold lined tickbox (). These will be delivered with the wheelchair if no other option is selected.
- For further information on different configurations and accessories, see the Prescriber Information (doc.nr. 95715-1) and each specific Info Sheet. User manuals and assembly instructions are also available at hdrehab.com.

1. CHASSIS

Info Sheet 95937-1 *Chassis*

| | Seat width | | Standard | | Reinforced |
|--|------------|--------------------------|----------|--------------------------|------------|
| HD Balance 24 Self propelled version Seat depth 43-52cm (with seat cushion) | 38 | <input type="checkbox"/> | 60210-6 | <input type="checkbox"/> | 60220-6 |
| | 42 | <input type="checkbox"/> | 60310-6 | <input type="checkbox"/> | 60320-6 |
| | 46 | <input type="checkbox"/> | 60410-6 | <input type="checkbox"/> | 60420-6 |
| | 50 | <input type="checkbox"/> | 60510-6 | <input type="checkbox"/> | 60520-6 |
| HD Balance 24 Self propelled version Extended seat depth 48-57cm (with seat cushion) | 38 | <input type="checkbox"/> | 60215-6 | <input type="checkbox"/> | 60225-6 |
| | 42 | <input type="checkbox"/> | 60315-6 | <input type="checkbox"/> | 60325-6 |
| | 46 | <input type="checkbox"/> | 60415-6 | <input type="checkbox"/> | 60425-6 |
| | 50 | <input type="checkbox"/> | 60515-6 | <input type="checkbox"/> | 60525-6 |
| HD Balance 16 Transport model Seat depth 43-52cm (with seat cushion) | 38 | <input type="checkbox"/> | 60250-6 | <input type="checkbox"/> | 60260-6 |
| | 42 | <input type="checkbox"/> | 60350-6 | <input type="checkbox"/> | 60360-6 |
| | 46 | <input type="checkbox"/> | 60450-6 | <input type="checkbox"/> | 60460-6 |
| | 50 | <input type="checkbox"/> | 60550-6 | <input type="checkbox"/> | 60560-6 |
| HD Balance 16 Transport model Extended seat depth 48-57cm (with seat cushion) | 38 | <input type="checkbox"/> | 60255-6 | <input type="checkbox"/> | 60265-6 |
| | 42 | <input type="checkbox"/> | 60355-6 | <input type="checkbox"/> | 60365-6 |
| | 46 | <input type="checkbox"/> | 60455-6 | <input type="checkbox"/> | 60465-6 |
| | 50 | <input type="checkbox"/> | 60555-6 | <input type="checkbox"/> | 60565-6 |



2. WHEELS & CASTORS

Info Sheet 95938-1 *Wheels*

| | <input type="checkbox"/> | HD Balance 24 | HD Balance 16 | |
|--|-------------------------------------|---------------|---------------|--|
| PUR (Polyurethane) drive wheels <i>Solid tyres with aluminium hand rims</i> | <input checked="" type="checkbox"/> | 67505-1 | 67551-1 | |
| Flexel drive wheels <i>Solid tyres with aluminium hand rims</i> | <input type="checkbox"/> | 67506-1 | 67552-1 | |
| Pneumatic drive wheels <i>With aluminium hand rims</i> | <input type="checkbox"/> | 67507-1 | 67553-1 | |

| | | | |
|--|-------------------------------------|-------|--|
| Castors 175 x 43mm | <input checked="" type="checkbox"/> | 67310 | |
| Castors 175 x 28mm | <input type="checkbox"/> | 67316 | |
| Castors 125 x 27mm, extended fork | <input type="checkbox"/> | 67318 | |
| Castors Flexel 200 x 50mm <i>Conflict with footrest in 90° position</i> | <input type="checkbox"/> | 67314 | |
| Castors 200 x 28mm <i>Conflict with footrest in 90° position</i> | <input type="checkbox"/> | 67317 | |

3. HAND RIMS (HD Balance 24)

Info Sheet 95938-1 *Wheels*

| | | | |
|---------------------------|-------------------------------------|-------|--|
| Aluminium | <input checked="" type="checkbox"/> | 15491 | |
| Stainless steel | <input type="checkbox"/> | 15492 | |
| Soft grip, plastic coated | <input type="checkbox"/> | 15943 | |

4. PUSH BAR / HANDLES

Info Sheet 95939-1 *Push bar / Handles / Caregiver brakes*

| | | | |
|---|-------------------------------------|------------------|--|
| Side unit, standard, left and right <i>With caregiver brakes</i> | <input checked="" type="checkbox"/> | 62330-1, 62330-2 | |
| Side unit, extended +12cm, left and right <i>With caregiver brakes. Longer wires for brakes, seat tilt, and back recline included.</i> | <input type="checkbox"/> | 62330-5, 62330-6 | |
| Side unit, angle adjustable, left and right <i>With caregiver brakes</i> | <input type="checkbox"/> | 62392-1, 62392-2 | |

| | <input type="checkbox"/> | sw38 | sw42 | sw46 | sw50 | |
|---|-------------------------------------|---------|---------|---------|---------|--|
| Push bar, standard | <input checked="" type="checkbox"/> | 62341-1 | 63341-1 | 64341-1 | 65341-1 | |
| Push bar, extra depth (+8cm) | <input type="checkbox"/> | 62342-1 | 63342-1 | 64342-1 | 65342-1 | |
| Push bar, extended & bent <i>Not in combination with Side units extended or angle adjustable</i> | <input type="checkbox"/> | 62343-1 | 63343-1 | 64343-1 | 65343-1 | |
| Separate operating handles <i>Not in combination with Side units angle adjustable</i> | <input type="checkbox"/> | 10295 | | | | |

5. SEAT CUSHION

Info Sheet 95940-1 *Seat cushions*


| | | sw38 | sw42 | sw46 | sw50 |
|---|------------------------------|---------|-------|-------|-------|
| Seat cushion standard, Grey Trevira | <input type="checkbox"/> | 68201 | 68301 | 68401 | 68501 |
| Seat cushion standard, Black Trevira | <input type="checkbox"/> | 68203 | 68303 | 68403 | 68503 |
| Seat cushion standard, Hygiene | <input type="checkbox"/> | 68204 | 68304 | 68404 | 68504 |
| Seat cushion long, +5cm, Grey Trevira | <input type="checkbox"/> | 68205 | 68305 | 68405 | 68505 |
| Seat cushion long, +5cm, Black Trevira | <input type="checkbox"/> | 68207 | 68307 | 68407 | 68507 |
| Seat cushion long, +5cm, Hygiene <i>(for chassis with extended seat depth)</i> | <input type="checkbox"/> | 68208 | 68308 | 68408 | 68508 |
| Additional cover, standard, Grey | <input type="checkbox"/> pcs | 49211 | 49311 | 49411 | 49511 |
| Additional cover, standard, Black | <input type="checkbox"/> pcs | 49213 | 49313 | 49413 | 49513 |
| Additional cover, standard, Hygiene | <input type="checkbox"/> pcs | 49214 | 49314 | 49414 | 49514 |
| Additional cover, long, Grey | <input type="checkbox"/> pcs | 49215 | 49315 | 49415 | 49515 |
| Additional cover, long, Black | <input type="checkbox"/> pcs | 49217 | 49317 | 49417 | 49517 |
| Additional cover, long, Hygiene | <input type="checkbox"/> pcs | 49218 | 49318 | 49418 | 49518 |
| Elastic seat cover | <input type="checkbox"/> pcs | 49205 | 49305 | 49405 | 49505 |
| Elastic seat cover, Incontinence | <input type="checkbox"/> pcs | 49206 | 49306 | 49406 | 49506 |
| Elastic seat cover +5cm | <input type="checkbox"/> pcs | 49207 | 49307 | 49407 | 49507 |
| Elastic seat cover +5cm, Incontinence | <input type="checkbox"/> pcs | 49208 | 49308 | 49408 | 49508 |
| <i>Elastic seat covers are fitted over existing cover. All elastic seat covers are Black Trevira.</i> | | | | | |
| Without seat cushion | <input type="checkbox"/> | 68200-8 | | | |



6. BACKREST SYSTEM

Info Sheet 95941-1 *Backrest systems*


























| | | sw38 | sw42 | sw46 | sw50 |
|-------------------------|--------------------------|-------|-------|-------|-------|
| Flexi-back standard | <input type="checkbox"/> | 62250 | 63250 | 64250 | 65250 |
| Flexi-back long (+10cm) | <input type="checkbox"/> | 62255 | 63255 | 64255 | 65255 |
| Solid back standard | <input type="checkbox"/> | 62260 | 63260 | 64260 | 65260 |
| Solid back long (+10cm) | <input type="checkbox"/> | 62265 | 63265 | 64265 | 65265 |





Note that increased or decreased width for the backrest system can be selected in section 15.

7. BACK CUSHION

Info Sheet 95941-1 *Backrest systems*

| | | sw38 | sw42 | sw46 | sw50 |
|---|---|---------|-------|-------|-------|
| Back cushion standard, Grey Trevira |  | 68281 | 68381 | 68481 | 68581 |
| Back cushion standard, Black Trevira |  | 68291 | 68391 | 68491 | 68591 |
| Back cushion standard, Dartex |  | 68282 | 68382 | 68482 | 68582 |
| <i>with side supporting wings</i> | | | | | |
| Back cushion flat, Grey Trevira |  | 68285 | 68385 | 68485 | 68585 |
| Back cushion flat, Black Trevira |  | 68295 | 68395 | 68495 | 68595 |
| Back cushion flat, Dartex |  | 68286 | 68386 | 68486 | 68586 |
| <i>without side supporting wings</i> | | | | | |
| Back cushion std long (+10cm), Grey |  | 68283 | 68383 | 68483 | 68583 |
| Back cushion std long (+10cm), Black |  | 68293 | 68393 | 68493 | 68593 |
| Back cushion std long (+10cm), Dartex |  | 68284 | 68384 | 68484 | 68584 |
| <i>with side supporting wings</i> | | | | | |
| Back cushion flat long (+10cm), Grey |  | 68287 | 68387 | 68487 | 68587 |
| Back cushion flat long (+10cm), Black |  | 68297 | 68397 | 68497 | 68597 |
| Back cushion flat long (+10cm), Dartex |  | 68288 | 68388 | 68488 | 68588 |
| <i>without side supporting wings</i> | | | | | |
| Additional cover standard, Grey Trevira |  pcs | 48221 | 48321 | 48421 | 48521 |
| Additional cover standard, Black Trevira |  pcs | 48231 | 48331 | 48431 | 48531 |
| Additional cover standard, Dartex |  pcs | 48222 | 48322 | 48422 | 48522 |
| <i>with side supporting wings</i> | | | | | |
| Additional cover flat, Grey Trevira |  pcs | 48225 | 48325 | 48425 | 48525 |
| Additional cover flat, Black Trevira |  pcs | 48235 | 48335 | 48435 | 48535 |
| Additional cover flat, Dartex |  pcs | 48226 | 48326 | 48426 | 48526 |
| <i>with side supporting wings</i> | | | | | |
| Additional cover std long, Grey Trevira |  pcs | 48223 | 48323 | 48423 | 48523 |
| Additional cover std long, Black Trevira |  pcs | 48233 | 48333 | 48433 | 48533 |
| Additional cover std long, Dartex |  pcs | 48224 | 48324 | 48424 | 48524 |
| <i>with side supporting wings</i> | | | | | |
| Additional cover flat long, Grey Trevira |  pcs | 48227 | 48327 | 48427 | 48527 |
| Additional cover flat long, Black Trevira |  pcs | 48237 | 48337 | 48437 | 48537 |
| Additional cover flat long, Dartex |  pcs | 48228 | 48328 | 48428 | 48528 |
| <i>without side supporting wings</i> | | | | | |
| Without back cushion |  | 68250-8 | | | |



8. LEG RESTS

Info Sheet 95933-1 Leg rests


| | Left | Right |
|--|---|---|
| Angle adjustable leg rest, standard | <input checked="" type="checkbox"/> 62725 | <input checked="" type="checkbox"/> 62720 |
| Angle adjustable leg rest, short | <input type="checkbox"/> 62721 | <input type="checkbox"/> 62722 |
| Angle adjustable leg rest, long | <input type="checkbox"/> 62723 | <input type="checkbox"/> 62724 |
| Disengage mechanism - knee joint <i>One sided adjustment of leg rest angle for one piece footrest / calf support / foot box</i> | <input type="checkbox"/> 62733-5 | <i>RH version on request</i> |
| Amputee leg rest <i>Angle adjustable, including calf support</i> | <input type="checkbox"/> sw38 - 62771 sw42 - 63771 sw46 - 64771 sw50 - 65771 | <input type="checkbox"/> sw38 - 62772 sw42 - 63772 sw46 - 64772 sw50 - 65772 |
| Fixed angle leg rest, 90°, standard | <input type="checkbox"/> 63758-1 | <input type="checkbox"/> 63758-2 |
| Fixed angle leg rest, 90°, short | <input type="checkbox"/> 63758-5 | <input type="checkbox"/> 63758-6 |
| Fixed angle leg rest, 90°, long | <input type="checkbox"/> 63758-3 | <input type="checkbox"/> 63758-4 |
| Fixed angle leg rest, 105°, standard | <input type="checkbox"/> 63759-1 | <input type="checkbox"/> 63759-2 |
| Fixed angle leg rest, 105°, short | <input type="checkbox"/> 63759-5 | <input type="checkbox"/> 63759-6 |
| Fixed angle leg rest, 105°, long | <input type="checkbox"/> 63759-3 | <input type="checkbox"/> 63759-4 |
| Without leg rest <i>Calf supports and Footrests then also excluded</i> | <input type="checkbox"/> 62701-8 | <input type="checkbox"/> 62702-8 |



9. CALF SUPPORTS












Info Sheet 95934-1 Calf supports

| | sw38 | sw42 | sw46 | sw50 |
|---|---|---|---|---|
| Calf support separate, left | <input checked="" type="checkbox"/> 62801 | <input checked="" type="checkbox"/> 63801 | <input checked="" type="checkbox"/> 64801 | <input checked="" type="checkbox"/> 65801 |
| Calf support separate, right | <input checked="" type="checkbox"/> 62802 | <input checked="" type="checkbox"/> 63802 | <input checked="" type="checkbox"/> 64802 | <input checked="" type="checkbox"/> 65802 |
| Calf support sep., +30mm extra depth, left | <input type="checkbox"/> 62803 | <input type="checkbox"/> 63803 | <input type="checkbox"/> 64803 | <input type="checkbox"/> 65803 |
| Calf support sep., +30mm extra depth, right | <input type="checkbox"/> 62804 | <input type="checkbox"/> 63804 | <input type="checkbox"/> 64804 | <input type="checkbox"/> 65804 |
| Calf support sep., +50mm extra depth, left | <input type="checkbox"/> 62805 | <input type="checkbox"/> 63805 | <input type="checkbox"/> 64805 | <input type="checkbox"/> 65805 |
| Calf support sep., +50mm extra depth, right | <input type="checkbox"/> 62806 | <input type="checkbox"/> 63806 | <input type="checkbox"/> 64806 | <input type="checkbox"/> 65806 |
| One piece calf support | <input type="checkbox"/> 62811 | <input type="checkbox"/> 63811 | <input type="checkbox"/> 64811 | <input type="checkbox"/> 65811 |
| One piece calf support, +30mm extra depth <i>Order disengage mechanism for knee joint separately</i> | <input type="checkbox"/> 62812 | <input type="checkbox"/> 63812 | <input type="checkbox"/> 64812 | <input type="checkbox"/> 65812 |
| Calf support panel | <input type="checkbox"/> 63891 | 65891 | | |
| Calf support strap | <input type="checkbox"/> 63892 | 63892 | | |
| Mounts for calf panel / calf strap, 2 pcs | <input type="checkbox"/> 63895 | | | |
| Without calf support left | <input type="checkbox"/> 63801-8 | | | |
| Without calf support right | <input type="checkbox"/> 63802-8 | | | |




















10. FOOTRESTS

Info Sheet 95935-1 Footrests

| | | sw38 | sw42 | sw46 | sw50 | |
|--|--------------------------|---------|---------|---------|---------|---|
| Separate footrest, standard, left | <input type="checkbox"/> | 62909 | 63909 | 64909 | 65909 |  |
| Separate footrest, standard, right | <input type="checkbox"/> | 62908 | 93908 | 64908 | 65908 | |
| Separate footrest, short, left | <input type="checkbox"/> | 62903 | 63903 | 64903 | 65903 |  |
| Separate footrest, short, right <i>For use with short leg rests</i> | <input type="checkbox"/> | 62904 | 63904 | 64904 | 65904 | |
| One-piece footrest, standard <i>Order disengage mechanism for knee joint separately</i> | <input type="checkbox"/> | 62951 | 63951 | 64951 | 65951 |  |
| One-piece footrest, short <i>For use with short leg rests</i> <i>Order disengage mechanism for knee joint separately</i> | <input type="checkbox"/> | 62953 | 63953 | 64953 | 65953 |  |
| Footbox, standard <i>Order disengage mechanism for knee joint separately</i> | <input type="checkbox"/> | 62981 | 63981 | 64981 | 65981 |  |
| Footbox, short <i>For use with short leg rests</i> <i>Order disengage mechanism for knee joint separately</i> | <input type="checkbox"/> | 62982 | 63982 | 64982 | 65982 |  |
| Padding, separate <i>Kit with 2 pieces, for separate footrests</i> | <input type="checkbox"/> | 62991-1 | | 64991-1 | |  |
| Padding, one-piece <i>For one-piece footrest or footbox</i> | <input type="checkbox"/> | 62993-1 | 63993-1 | 64993-1 | 65993-1 |  |
| Footrest side support, left | <input type="checkbox"/> | 65987 | | | |  |
| Footrest side support, right <i>For use with standard length footrests</i> | <input type="checkbox"/> | 65986 | | | | |
| Mounting of footrest side support | <input type="checkbox"/> | 65987-8 | | | | |
| Heel strap, kit with 2 pieces <i>For use with separate footrests</i> | <input type="checkbox"/> | 62995-1 | | | |  |
| Mounting of heel straps | <input type="checkbox"/> | 62995-8 | | | | |
| Without footrest left | <input type="checkbox"/> | 63909-8 | | | |  |
| Without footrest right | <input type="checkbox"/> | 63908-8 | | | | |

11. ARMRESTS

Info Sheet 95936-1 Armrests

| | Left | Right | |
|--|----------------------------------|----------------------------------|---|
| Armrest standard | <input type="checkbox"/> 62401 | <input type="checkbox"/> 62402 |  |
| Armrest wide | <input type="checkbox"/> 63401 | <input type="checkbox"/> 63402 |  |
| Armrest low | <input type="checkbox"/> 62405 | <input type="checkbox"/> 62406 |  |
| Armrest wide low | <input type="checkbox"/> 63405 | <input type="checkbox"/> 63406 |  |
| Armrest standard with wedge side support 7° | <input type="checkbox"/> 62411 | <input type="checkbox"/> 62412 |  |
| Armrest standard with wedge side support 4° | <input type="checkbox"/> 62413 | <input type="checkbox"/> 62414 |  |
| Armrest wide with wedge side support 7° | <input type="checkbox"/> 63411 | <input type="checkbox"/> 63412 |  |
| Armrest wide with wedge side support 4° | <input type="checkbox"/> 63413 | <input type="checkbox"/> 63414 |  |
| Hemiplegia armrest, high | <input type="checkbox"/> 63491 | <input type="checkbox"/> 63492 |  |
| Hemiplegia armrest, low | <input type="checkbox"/> 63493 | <input type="checkbox"/> 63494 |  |
| Side support standard with post, w/o arm pad | <input type="checkbox"/> 62451-1 | <input type="checkbox"/> 62451-2 |  |
| Post with cap, without arm pad | <input type="checkbox"/> 62453-1 | <input type="checkbox"/> 62453-1 |  |
| Wedge side support 7° | <input type="checkbox"/> 63463 | <input type="checkbox"/> 63464 |  |
| Wedge side support 4° | <input type="checkbox"/> 63461 | <input type="checkbox"/> 63462 |  |
| Coaming for armrest, standard | <input type="checkbox"/> 62471 | <input type="checkbox"/> 62472 |  |
| Coaming for armrest, wide | <input type="checkbox"/> 63471 | <input type="checkbox"/> 63472 |  |
| Mounting of armrest coaming | <input type="checkbox"/> 62471-8 | |  |
| Armrest locks, pair <i>Locks the armrests to the seat frame</i> | <input type="checkbox"/> 62459 | | |
| Mounting of armrest locks | <input type="checkbox"/> 62459-8 | | |
| Without armrest | <input type="checkbox"/> 62400-8 | <input type="checkbox"/> 62400-8 | |

ALTERNATE CONFIGURATIONS

Alternate configurations of HD Balance wheelchairs must be evaluated for risk. Evaluate if there is an increased risk of tipping and/or pinching for the particular user with the desired modification. If there is increased risk, make an informed decision as to whether the alternate configuration is still warranted. **Inform everyone involved about the risks.**

12. BALANCE POINT

Info Sheet 95942-1 Balance Point

| | | | |
|--|--------------------------|---------|---|
| Forward mounted seat unit | <input type="checkbox"/> | 67056-8 | See 95758-1 <i>Balancing modes</i> , 95759-1 <i>Reference table</i> |
| Back mounted seat unit | <input type="checkbox"/> | 67057-8 | <i>balancing configurations</i> , and 95707-1 <i>Risk information</i> |
| Drive wheels rear position (5 cm back) | <input type="checkbox"/> | 67551-8 | |

13. SEAT HEIGHT AND SEAT DEPTH

Info Sheet 95943-1 Seat height, 95937-1 Chassis

| | | | |
|-----------------------------|--------------------------|---------|--|
| Increased seat height +5 cm | <input type="checkbox"/> | 67052-8 | See 95758-1 <i>Balancing modes</i> , 95759-1 <i>Reference table</i> <i>balancing configurations</i> , and 95707-1 <i>Risk information</i> |
| Lowered seat height -3 cm | <input type="checkbox"/> | 67053-8 | |

| | | | |
|---|-----------|--------------------------|---------|
| 20" wheels for HD Balance 24 <i>Lowers seat height approx. 4 cm</i> | PUR | <input type="checkbox"/> | 67515-2 |
| | Flexel | <input type="checkbox"/> | 67516-2 |
| | Pneumatic | <input type="checkbox"/> | 67518-2 |
| 125x27 castors for HD Balance 24 <i>Only used in combination with 20" wheels</i> | | <input type="checkbox"/> | 67315 |



| | | | |
|---|--------------------------|---------|-------------------------------------|
| Back spacers, kit with 4 pieces <i>Can be mounted on one side if desired</i> | <input type="checkbox"/> | 69791-5 | See 95707-1 <i>Risk information</i> |
| Mounting of back spacers | <input type="checkbox"/> | 69791-8 | |



14. SEAT TILT AND BACK RECLINE

Info Sheet 95944-1 Seat tilt and Back recline

See 95707-1 *Risk information*. Standard configuration has 0° - 20° seat tilt and 90° - 120° back recline.

| | | | |
|--|--------------------------|---|--------------------------|
| Seat tilt -5° - 20° <i>Allows 5° forward tilt</i> | <input type="checkbox"/> | 67820-8 | |
| Increased seat tilt 0° - 30° <i>Note that drive wheels should be in rear position</i> | <input type="checkbox"/> | 67823-8 | |
| Increased seat tilt -5° - 30° <i>Allows 5° forward tilt</i> <i>Note that drive wheels should be in rear position</i> | <input type="checkbox"/> | 67824-8 | |
| Increased back recline 90° - 135° <i>Note that drive wheels should be in rear position</i> | <input type="checkbox"/> | 67705-8 | |
| Lockable back recline angle | <input type="checkbox"/> | 62369 | |
| Fixed back recline angle <i>Control and gas spring to reset are not included</i> | <input type="checkbox"/> | Recline angle on delivery (90° - 105° possible): | <input type="checkbox"/> |

15. INCREASED OR DECREASED BACK WIDTH

Info Sheet 95945-1 *Back System - increased or decreased width*

| | sw38 | sw42 | sw46 | sw50 |
|---|------------|---------|------------|------------|
| Wider back system +4 cm, incl. wider back cushion <input type="checkbox"/> | 63298-8 | 64298-8 | 65298-8 | On request |
| Wider back system +8 cm, incl. wider back cushion <input type="checkbox"/> <i>Causes certain limitations for the backrest height setting and armrest depth</i> | 64299-8 | 65299-8 | On request | On request |
| Narrower back system -4 cm, incl. narrower cushion <input type="checkbox"/> | On request | 62297-8 | 63297-8 | 64297-8 |

16. CONFIGURATION FOR SEVERE LOAD

Info Sheet 95946-1 *Configurations for severe load from user*

Shock absorbing back
Available for both adjustable and fixed back recline angle

17. BRAKES, SPECIAL VERSIONS

Info Sheet 95947-1 *Brakes, Special versions*

User brakes for HD Balance 16 67130-1
 User brakes with extended lever 67115-5
For HD Balance 24

18. OUTER DIMENSIONS

Info Sheet 95938-1 *Wheels, 95719-1 Measures*

Castor pin narrow 69991-1
Reduces total width (in combination with narrow position for hand rims)
 Narrow position for hand rims 15499-8
Reduces total width, not possible for user to drive

ACCESSORIES

19. HEADREST / NECK REST

















Info Sheet 95932-1 Headrest / Neck rest

| | Grey Trevira | Black Trevira | Dartex |
|--|--------------------------------------|--------------------------------------|--------------------------------------|
| Headrest Basic | <input type="checkbox"/> 69108-1 | <input type="checkbox"/> 69108-3 | <input type="checkbox"/> 69108-2 |
| Headrest Side, sideways adj. | <input type="checkbox"/> 69101-1 | <input type="checkbox"/> 69101-3 | <input type="checkbox"/> 69101-2 |
| Headrest Side+ <i>With mounts for side supports Side supports ordered separately</i> | <input type="checkbox"/> 69104-1 | <input type="checkbox"/> 69104-3 | <input type="checkbox"/> 69104-2 |
| Side support for headrest, Vario <i>Adjustable angle</i> | <input type="checkbox"/> pcs 69182-1 | <input type="checkbox"/> pcs 69182-3 | <input type="checkbox"/> pcs 69182-2 |
| Side support for headrest, Fix <i>Fixed angle</i> | <input type="checkbox"/> pcs 69181-1 | <input type="checkbox"/> pcs 69181-3 | <input type="checkbox"/> pcs 69181-2 |
| Headrest cover <i>Black Trevira</i> | | <input type="checkbox"/> pcs 49911 | |
| Neck rest Side | <input type="checkbox"/> 69111-1 | <input type="checkbox"/> 69111-3 | <input type="checkbox"/> 69111-2 |
| Neck rest Side+ <i>With mounts for side supports Side supports ordered separately</i> | <input type="checkbox"/> 69114-1 | <input type="checkbox"/> 69114-3 | <input type="checkbox"/> 69114-2 |
| Side support for neck rest, Vario <i>Adjustable angle</i> | <input type="checkbox"/> pcs 69191-1 | <input type="checkbox"/> pcs 69191-3 | <input type="checkbox"/> pcs 69191-2 |
| Extra cover, Headrest Basic | <input type="checkbox"/> pcs 49922 | <input type="checkbox"/> pcs 49952 | <input type="checkbox"/> pcs 49924 |
| Extra cover, Headrest Side | <input type="checkbox"/> pcs 49921 | <input type="checkbox"/> pcs 49951 | <input type="checkbox"/> pcs 49923 |
| Extra cover, Neck rest Side | <input type="checkbox"/> pcs 49926 | <input type="checkbox"/> pcs 49953 | <input type="checkbox"/> pcs 49927 |
| Extra cover, Side support <i>For Headrest and Neck rest</i> | <input type="checkbox"/> pcs 49931 | <input type="checkbox"/> pcs 49941 | <input type="checkbox"/> pcs 49932 |
| Special variant: mounting bar extended +10cm | | <input type="checkbox"/> 12780-1 | |



20. THORACIC SUPPORT & SHOULDER SUPPORT


Info Sheet 95931-1 *Thoracic support*, 95954-1 *Shoulder support*

| | | Grey Trevira | Black Trevira | Dartex | | | | |
|--|--------------------------|--|----------------------------|--|----------------------------|--|---|---|
| Thoracic supp. Standard, narrow pos. | <input type="checkbox"/> | pcs 69301-1 | <input type="checkbox"/> | pcs 69301-3 | <input type="checkbox"/> | pcs 69301-2 |  | |
| Thoracic supp. Standard, wide pos. | <input type="checkbox"/> | pcs 69302-1 | <input type="checkbox"/> | pcs 69302-3 | <input type="checkbox"/> | pcs 69302-2 |  | |
| Thoracic support Adjustable Short | <input type="checkbox"/> | pcs 69311-1 | <input type="checkbox"/> | pcs 69311-3 | <input type="checkbox"/> | pcs 69311-2 |  | |
| Thoracic support Adjustable Long | <input type="checkbox"/> | pcs 69312-1 | <input type="checkbox"/> | pcs 69312-3 | <input type="checkbox"/> | pcs 69312-2 |  | |
| Thoracic support Adjustable Mini | | | <input type="checkbox"/> | pcs 69332-3 | | |  | |
| Thoracic support Standard, narrow position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69305-1 pcs 69306-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69305-3 pcs 69306-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69305-2 pcs 69306-2 |  |
| Thoracic support Standard, wide position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69307-1 pcs 69308-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69307-3 pcs 69308-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69307-2 pcs 69308-2 |  |
| Thoracic supp. Adj. Short, narrow position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69315-1 pcs 69316-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69315-3 pcs 69316-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69315-2 pcs 69316-2 |  |
| Thoracic supp. Adj. Short, wide position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69317-1 pcs 69318-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69317-3 pcs 69318-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69317-2 pcs 69318-2 |  |
| Thoracic supp. Adj. Long, narrow position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69325-1 pcs 69326-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69325-3 pcs 69326-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69325-2 pcs 69326-2 |  |
| Thoracic supp. Adj. Long, wide position, with Swing | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69327-1 pcs 69328-1 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69327-3 pcs 69328-3 | <input type="checkbox"/> <input type="checkbox"/> | pcs 69327-2 pcs 69328-2 |  |
| Thoracic supp. Adj. Mini, narrow position, with Swing | | | left right | <input type="checkbox"/> <input type="checkbox"/> | pcs 69323-1 pcs 69323-2 | | |  |
| Thoracic support mount | <input type="checkbox"/> | pcs 69351-1 | | | | | |  |
| Mounting of thoracic support mount | <input type="checkbox"/> | 69351-8 | | | | | |  |
| Shoulder support, Swing, complete with mount | | | left right | <input type="checkbox"/> <input type="checkbox"/> | 69307-3 69308-3 | | |  |
| Mounting of shoulder support | <input type="checkbox"/> | 69201-8 | | | | | |  |

21. POMMELS

Info Sheet 95948-1 *Pommels*

| | | |
|--------------------------|--------------------------|---------|
| Leg divider, small | <input type="checkbox"/> | 69605-2 |
| Leg divider, large | <input type="checkbox"/> | 69606-2 |
| Knee divider, small | <input type="checkbox"/> | 69615-2 |
| Knee divider, large | <input type="checkbox"/> | 69616-2 |
| Pommel mount | <input type="checkbox"/> | 69620-1 |
| Mounting of pommel mount | <input type="checkbox"/> | 69620-8 |



22. TRAY TABLE

Info Sheet 95930-1 *Tray tables*

| | | | | | |
|------------------------------------|-------------------------------|-------------------|--------------------------------|-------------------|---------------|
| Tray table, PC (smoke coloured) | <input type="checkbox"/> | sw38 69521 | sw42 69531 | sw46 69541 | sw50 69551 |
| Tray table, PC (smoke), with pads | <input type="checkbox"/> | 69522 | 69532 | 69542 | 69552 |
| Tray table, Wood (oak) | <input type="checkbox"/> | 69525 | 69535 | 69545 | 69555 |
| Tray table, Wood (oak), with pads | <input type="checkbox"/> | 69526 | 69536 | 69546 | 69556 |
| Support table, padded | <input type="checkbox"/> | ----- 69596 ----- | | ----- 69598 ----- | |
| Mug tray (oak) | left <input type="checkbox"/> | 69581 | right <input type="checkbox"/> | 69582 | |
| Tray table mounts, pair | <input type="checkbox"/> | 69515-1 | | | |
| Tray table mounts w/ locking, pair | <input type="checkbox"/> | 69515-5 | | | |
| Mug tray mount w/ locking | left <input type="checkbox"/> | 69585-1 | right <input type="checkbox"/> | 69585-2 | |
| Mounting of tray mounts | <input type="checkbox"/> | 69515-8 | | | |



23. POSITIONING BELT & MOUNTS

Info Sheet 95906-1 *Pos. belt*, 95949-1 *Belt mounts*, 95959-1 *Harness mounts*

2-point Positioning belt
Order 1 set of belt mounts separately

13545



4-point Positioning belt
Order 2 sets of belt mounts separately

13546



Belt mounts, set with 2 pieces

sets 69710-1



Mounting of belt mounts

69710-8

Harness mounts for Flexi-back, pair
Including strap holders for lower straps

69720-5



Harness mounts for Solid back, pair
Including strap holders for lower straps

69730-5



Mounting of harness mounts

69720-8 (Flexi)
69730-8 (Solid)


24. SPOKE PROTECTION & WHEEL COVERS

Info Sheet 95938-1 *Wheels*

Spoke protection for HD Balance 24, pair (24" - 15471, 20" - 15472)

Mounting of spoke protection

15475-8



Wheel covers for HD Balance 16, pair 15690

Mounting of wheel covers

15690-8



25. ELECTRICAL FUNCTIONS

Info Sheet 95951-1 *Electrical functions*

Electric seat tilt and back recline
Only in combination with standard seat tilt (0-20°) and back recline (90-120°)









69801



Electric seat tilt
Only in combination with standard seat tilt (0-20°)

69802

26. OTHER ACCESSORIES

| | | | |
|---|--------------------------------|-----------|--|
| Side pads 1 cm, pair <i>Reduces seat width 1 cm per side. Info Sheet 95937-1.</i> | <input type="checkbox"/> | 69680 |  |
| Spacers for side supports 1 cm, pair <i>Reduces seat width 1 cm per side. Info Sheet 95937-1.</i> | <input type="checkbox"/> | 62458-1 |  |
| Back wedge <i>Pushes back cushion forward on one side. Info Sheet 95958-1.</i> | <input type="checkbox"/> | 69698-1 |  |
| Profile stop <i>Assists in mounting supports (head/neck/thoracic) in correct position. Info Sheet 95932.</i> | <input type="checkbox"/> | pcs 12799 |  |
| Bag hook <i>Info Sheet 95957-1.</i> | <input type="checkbox"/> | pcs 13605 |  |
| Mounting of bag hook | <input type="checkbox"/> | 13605-8 | |
| IV pole <i>Info Sheet 95956-1 and Risk Information 95707-1.</i> | <input type="checkbox"/> | 69760 |  |
| IV pole mount | left <input type="checkbox"/> | 69755-1 |  |
| | right <input type="checkbox"/> | 69755-2 | |
| Mounting of IV pole mount | <input type="checkbox"/> | 69750-8 | |
| Oxygen tank holder 2 litre <i>Info Sheet 95955-1 and Risk Information 95707-1.</i> | <input type="checkbox"/> | 59605 |  |
| Mounting of oxygen tank holder | <input type="checkbox"/> | 56904-8 | |

27. NOTES

