

Date	<input type="text"/>	Order number	<input type="text"/>
Customer	<input type="text"/>	Contact person	<input type="text"/>
Customer tel.	<input type="text"/>	Contact person tel.	<input type="text"/>
Customer email	<input type="text"/>	Contact person email	<input type="text"/>
User	<input type="text"/>	Contact at HD Rehab	<input type="text"/>

- This Prescription Form is used when ordering a **single** chassis (with accessories) in a particular configuration.
- All parts that are included in standard configuration are marked with a bold lined tickbox (). These will be delivered with the wheelchair if no other option is selected.
- For further information on different configurations and accessories, see the Prescriber Information (doc.nr. 95715-1) and each specific Info Sheet. User manuals and assembly instructions are also available at hdrehab.com.

1. CHASSIS

Info Sheet 95937-1 Chassis

	Seat width		Standard
Wheelchair chassis 24 Self propelled version	38	<input type="checkbox"/>	60230-6
	42	<input type="checkbox"/>	60330-6
	46	<input type="checkbox"/>	60430-6
	50	<input type="checkbox"/>	60530-6
Wheelchair chassis 24 Self propelled version Extended seat frame +5 cm	38	<input type="checkbox"/>	60235-6
	42	<input type="checkbox"/>	60335-6
	46	<input type="checkbox"/>	60435-6
	50	<input type="checkbox"/>	60535-6
Wheelchair chassis 16 Transport model Seat depth 43-52cm (with seat cushion)	38	<input type="checkbox"/>	60270-6
	42	<input type="checkbox"/>	60370-6
	46	<input type="checkbox"/>	60470-6
	50	<input type="checkbox"/>	60570-6
Wheelchair chassis 16 Transport model Extended seat frame +5 cm	38	<input type="checkbox"/>	60275-6
	42	<input type="checkbox"/>	60375-6
	46	<input type="checkbox"/>	60475-6
	50	<input type="checkbox"/>	60575-6



2. WHEELS & CASTORS

Info Sheet 95938-1 *Wheels*

		model 24	model 16
PUR (Polyurethane) drive wheels <i>Solid tyres with aluminium hand rims</i>	<input type="checkbox"/>	67505-1	67551-1
Flexel drive wheels <i>Solid tyres with aluminium hand rims</i>	<input type="checkbox"/>	67506-1	67552-1
Pneumatic drive wheels <i>With aluminium hand rims</i>	<input type="checkbox"/>	67507-1	67553-1




Castors 175 x 43mm	<input checked="" type="checkbox"/>	67310
Castors 175 x 28mm	<input type="checkbox"/>	67316
Castors 125 x 27mm, extended fork	<input type="checkbox"/>	67318
Castors Flexel 200 x 50mm <i>Conflict with footrest in 90° position</i>	<input type="checkbox"/>	67314
Castors 200 x 28mm <i>Conflict with footrest in 90° position</i>	<input type="checkbox"/>	67317



3. HAND RIMS (model 24)

Info Sheet 95938-1 *Wheels*


Aluminium	<input checked="" type="checkbox"/>	15491
Stainless steel	<input type="checkbox"/>	15492
Soft grip, plastic coated	<input type="checkbox"/>	15943



4. PUSH BAR / HANDLES

Info Sheet 95939-1 *Push bar / Handles / Caregiver brakes*

Side unit, standard, left and right <i>With caregiver brakes</i>	<input checked="" type="checkbox"/>	62330-1, 62330-2
Side unit, extended +12cm, left and right <i>With caregiver brakes. Longer wires for brakes, seat tilt, and back recline included.</i>	<input type="checkbox"/>	62330-5, 62330-6
Side unit, angle adjustable, left and right <i>With caregiver brakes</i>	<input type="checkbox"/>	62392-1, 62392-2



		sw38	sw42	sw46	sw50
Push bar, standard	<input checked="" type="checkbox"/>	62341-1	63341-1	64341-1	65341-1
Push bar, extra depth (+8cm)	<input type="checkbox"/>	62342-1	63342-1	64342-1	65342-1
Push bar, extended & bent <i>Not in combination with Side units extended or angle adjustable</i>	<input type="checkbox"/>	62343-1	63343-1	64343-1	65343-1
Separate operating handles <i>Not in combination with Side units angle adjustable</i>	<input type="checkbox"/>	10295			



5. LEG RESTS

Info Sheet 95933-1 Leg rests

	Left	Right
Angle adjustable leg rest, standard	<input checked="" type="checkbox"/> 62725	<input checked="" type="checkbox"/> 62720
Angle adjustable leg rest, short	<input type="checkbox"/> 62721	<input type="checkbox"/> 62722
Angle adjustable leg rest, long	<input type="checkbox"/> 62723	<input type="checkbox"/> 62724
Disengage mechanism - knee joint <i>One sided adjustment of leg rest angle for one piece footrest / calf support / foot box</i>	<input type="checkbox"/> 62733-5	<i>RH version on request</i>
Amputee leg rest <i>Angle adjustable, including calf support</i>	<input type="checkbox"/> sw38 - 62771 sw42 - 63771 sw46 - 64771 sw50 - 65771	<input type="checkbox"/> sw38 - 62772 sw42 - 63772 sw46 - 64772 sw50 - 65772
Fixed angle leg rest, 90°, standard	<input type="checkbox"/> 63758-1	<input type="checkbox"/> 63758-2
Fixed angle leg rest, 90°, short	<input type="checkbox"/> 63758-5	<input type="checkbox"/> 63758-6
Fixed angle leg rest, 90°, long	<input type="checkbox"/> 63758-3	<input type="checkbox"/> 63758-4
Fixed angle leg rest, 105°, standard	<input type="checkbox"/> 63759-1	<input type="checkbox"/> 63759-2
Fixed angle leg rest, 105°, short	<input type="checkbox"/> 63759-5	<input type="checkbox"/> 63759-6
Fixed angle leg rest, 105°, long	<input type="checkbox"/> 63759-3	<input type="checkbox"/> 63759-4
Without leg rest <i>Calf supports and Footrests then also excluded</i>	<input type="checkbox"/> 62701-8	<input type="checkbox"/> 62702-8



6. CALF SUPPORTS












Info Sheet 96934-1 Calf supports

	sw38	sw42	sw46	sw50
Calf support separate, left	<input checked="" type="checkbox"/> 62801	<input checked="" type="checkbox"/> 63801	<input checked="" type="checkbox"/> 64801	<input checked="" type="checkbox"/> 65801
Calf support separate, right	<input checked="" type="checkbox"/> 62802	<input checked="" type="checkbox"/> 63802	<input checked="" type="checkbox"/> 64802	<input checked="" type="checkbox"/> 65802
Calf support sep., +30mm extra depth, left	<input type="checkbox"/> 62803	<input type="checkbox"/> 63803	<input type="checkbox"/> 64803	<input type="checkbox"/> 65803
Calf support sep., +30mm extra depth, right	<input type="checkbox"/> 62804	<input type="checkbox"/> 63804	<input type="checkbox"/> 64804	<input type="checkbox"/> 65804
Calf support sep., +50mm extra depth, left	<input type="checkbox"/> 62805	<input type="checkbox"/> 63805	<input type="checkbox"/> 64805	<input type="checkbox"/> 65805
Calf support sep., +50mm extra depth, right	<input type="checkbox"/> 62806	<input type="checkbox"/> 63806	<input type="checkbox"/> 64806	<input type="checkbox"/> 65806
One piece calf support	<input type="checkbox"/> 62811	<input type="checkbox"/> 63811	<input type="checkbox"/> 64811	<input type="checkbox"/> 65811
One piece calf support, +30mm extra depth <i>Order disengage mechanism for knee joint separately</i>	<input type="checkbox"/> 62812	<input type="checkbox"/> 63812	<input type="checkbox"/> 64812	<input type="checkbox"/> 65812
Calf support panel	<input type="checkbox"/> 63891	65891		
Calf support strap	<input type="checkbox"/> 63892	63892		
Mounts for calf panel / calf strap, 2 pcs	<input type="checkbox"/> 63895			
Without calf support left	<input type="checkbox"/> 63801-8			
Without calf support right	<input type="checkbox"/> 63802-8			













7. FOOTRESTS

Info Sheet 96935-1 Footrests

		sw38	sw42	sw46	sw50	
Separate footrest, standard, left	<input type="checkbox"/>	62909	63909	64909	65909	
Separate footrest, standard, right	<input type="checkbox"/>	62908	93908	64908	65908	
Separate footrest, short, left	<input type="checkbox"/>	62903	63903	64903	65903	
Separate footrest, short, right <i>For use with short leg rests</i>	<input type="checkbox"/>	62904	63904	64904	65904	
One-piece footrest, standard <i>Order disengage mechanism for knee joint separately</i>	<input type="checkbox"/>	62951	63951	64951	65951	
One-piece footrest, short <i>For use with short leg rests</i> <i>Order disengage mechanism for knee joint separately</i>	<input type="checkbox"/>	62953	63953	64953	65953	
Footbox, standard <i>Order disengage mechanism for knee joint separately</i>	<input type="checkbox"/>	62981	63981	64981	65981	
Footbox, short <i>For use with short leg rests</i> <i>Order disengage mechanism for knee joint separately</i>	<input type="checkbox"/>	62982	63982	64982	65982	
Padding, separate <i>Kit with 2 pieces, for separate footrests</i>	<input type="checkbox"/>	62991-1		64991-1		
Padding, one-piece <i>For one-piece footrest or footbox</i>	<input type="checkbox"/>	62993-1	63993-1	64993-1	65993-1	
Footrest side support, left	<input type="checkbox"/>	65987				
Footrest side support, right <i>For use with standard length footrests</i>	<input type="checkbox"/>	65986				
Mounting of footrest side support	<input type="checkbox"/>	65987-8				
Heel strap, kit with 2 pieces <i>For use with separate footrests</i>	<input type="checkbox"/>	62995-1				
Mounting of heel straps	<input type="checkbox"/>	62995-8				
Without footrest left	<input type="checkbox"/>	63909-8				
Without footrest right	<input type="checkbox"/>	63908-8				

8. ARMRESTS

Info Sheet 95936-1 Armrests

	Left	Right	
Armrest standard	<input type="checkbox"/> 62401	<input type="checkbox"/> 62402	
Armrest wide	<input type="checkbox"/> 63401	<input type="checkbox"/> 63402	
Armrest low	<input type="checkbox"/> 62405	<input type="checkbox"/> 62406	
Armrest wide low	<input type="checkbox"/> 63405	<input type="checkbox"/> 63406	
Hemiplegia armrest, high	<input type="checkbox"/> 63491	<input type="checkbox"/> 63492	
Hemiplegia armrest, low	<input type="checkbox"/> 63493	<input type="checkbox"/> 63494	
Side support standard with post, w/o arm pad	<input type="checkbox"/> 62451-1	<input type="checkbox"/> 62451-2	
Coaming for armrest, standard	<input type="checkbox"/> 62471	<input type="checkbox"/> 62472	
Coaming for armrest, wide	<input type="checkbox"/> 63471	<input type="checkbox"/> 63472	
Mounting of armrest coaming	<input type="checkbox"/> 62471-8		
Armrest locks, pair <i>Locks the armrests to the seat frame</i>	<input type="checkbox"/> 62459		
Mounting of armrest locks	<input type="checkbox"/> 62459-8		
Without armrest	<input type="checkbox"/> 62400-8	<input type="checkbox"/> 62400-8	

ALTERNATE CONFIGURATIONS

Alternate configurations of HD Balance chassis must be evaluated for risk. Evaluate if there is an increased risk of tipping and/or pinching for the particular user with the desired modification. If there is increased risk, make an informed decision as to whether the alternate configuration is still warranted. **Inform everyone involved about the risks.**


9. BALANCE POINT

Info Sheet 96942-1 Balance Point

- | | | | |
|--|--------------------------|---------|--|
| Forward mounted seat unit | <input type="checkbox"/> | 67056-8 | See 95758-1 <i>Balancing modes</i> , 95759-1 <i>Reference table balancing configurations</i> , and 95707-1 <i>Risk information</i> |
| Back mounted seat unit | <input type="checkbox"/> | 67057-8 | |
| Drive wheels rear position (5 cm back) | <input type="checkbox"/> | 67551-8 | |

10. SEAT HEIGHT AND SEAT DEPTH

Info Sheet 95943-1 *Seat height*, 95937-1 *Chassis*

- | | | | |
|---|---------------------------------|---------|--|
| Increased seat height +5 cm | <input type="checkbox"/> | 67052-8 | See 95758-1 <i>Balancing modes</i> , 95759-1 <i>Reference table balancing configurations</i> , and 95707-1 <i>Risk information</i> |
| Lowered seat height -3 cm | <input type="checkbox"/> | 67053-8 | |
| 20" wheels for model 24 | <input type="checkbox"/> | 67515-2 |  |
| <i>Lowers seat height approx. 4 cm</i> | | | |
| | PUR <input type="checkbox"/> | 67516-2 | |
| | Flexel <input type="checkbox"/> | 67518-2 | |
| 125x27 castors for model 24 | <input type="checkbox"/> | 67315 | |
| <i>Only used in combination with 20" wheels</i> | | | |

11. SEAT TILT AND BACK RECLINE

Info Sheet 95944-1 *Seat tilt and Back recline*

See 95707-1 *Risk information*. Standard configuration has fixed back recline and 0° - 20° seat tilt.

- | | | | |
|---|-------------------------------------|-----------------------------|--------------------------|
| Fixed back recline angle (standard) | <input checked="" type="checkbox"/> | Recline angle on delivery | <input type="checkbox"/> |
| <i>Delivered with 95° if no selection is made</i> | | <i>(90°-105° possible):</i> | |
| Adjustable back recline | <input type="checkbox"/> | 67701-8 | |
| <i>Gas spring 90°-120°</i> | | | |
| Seat tilt -5° - 20° | <input type="checkbox"/> | 67820-8 | |
| <i>Allows 5° forward tilt</i> | | | |
| Increased seat tilt 0° - 30° | <input type="checkbox"/> | 67823-8 | |
| Note: <i>drive wheels should be in rear position</i> | | | |
| Increased seat tilt -5° - 30° | <input type="checkbox"/> | 67824-8 | |
| <i>Allows 5° forward tilt</i> | | | |
| Note: <i>drive wheels should be in rear position</i> | | | |

12. OUTER DIMENSIONS

Info Sheet 95938-1 *Wheels*, 95719-1 *Measures*

- | | | |
|--|--------------------------|---------|
| Castor pin narrow | <input type="checkbox"/> | 69991-1 |
| <i>Reduces total width (in combination with narrow position for hand rims)</i> | | |
| Narrow position for hand rims | <input type="checkbox"/> | 15499-8 |
| <i>Reduces total width, not possible for user to drive</i> | | |

13. CONFIGURATION FOR SEVERE LOAD

Info Sheet 96946-1 Configurations for severe load from user

Shock absorbing back
Available for both adjustable and fixed back recline angle

14. BRAKES, SPECIAL VERSIONS

Info Sheet 96947-1 Brakes, Special versions

User brakes for model16 67130-1
 User brakes with extended lever 67115-5
For model 24

ACCESSORIES

15. HEADREST / NECK REST

Info Sheet 95932-1 Headrest / Neck rest

	Grey Trevira	Black Trevira	Hygiene material
Headrest Basic	<input type="checkbox"/> 69108-1	<input type="checkbox"/> 69108-3	<input type="checkbox"/> 69108-2
Headrest Side, sideways adj.	<input type="checkbox"/> 69101-1	<input type="checkbox"/> 69101-3	<input type="checkbox"/> 69101-2
Headrest Side+ <i>With mounts for side supports Side supports ordered separately</i>	<input type="checkbox"/> 69104-1	<input type="checkbox"/> 69104-3	<input type="checkbox"/> 69104-2
Side support for headrest, Vario <i>Adjustable angle</i>	<input type="checkbox"/> pcs 69182-1	<input type="checkbox"/> pcs 69182-3	<input type="checkbox"/> pcs 69182-2
Side support for headrest, Fix <i>Fixed angle</i>	<input type="checkbox"/> pcs 69181-1	<input type="checkbox"/> pcs 69181-3	<input type="checkbox"/> pcs 69181-2
Headrest cover <i>Black Trevira</i>		<input type="checkbox"/> pcs 49911	
Neck rest Side	<input type="checkbox"/> 69111-1	<input type="checkbox"/> 69111-3	<input type="checkbox"/> 69111-2
Neck rest Side+ <i>With mounts for side supports Side supports ordered separately</i>	<input type="checkbox"/> 69114-1	<input type="checkbox"/> 69114-3	<input type="checkbox"/> 69114-2
Side support for neck rest, Vario <i>Adjustable angle</i>	<input type="checkbox"/> pcs 69191-1	<input type="checkbox"/> pcs 69191-3	<input type="checkbox"/> pcs 69191-2
Extra cover, Headrest Basic	<input type="checkbox"/> pcs 49922	<input type="checkbox"/> pcs 49952	<input type="checkbox"/> pcs 49924
Extra cover, Headrest Side	<input type="checkbox"/> pcs 49921	<input type="checkbox"/> pcs 49951	<input type="checkbox"/> pcs 49923
Extra cover, Neck rest Side	<input type="checkbox"/> pcs 49926	<input type="checkbox"/> pcs 49953	<input type="checkbox"/> pcs 49927
Extra cover, Side support <i>For Headrest and Neck rest</i>	<input type="checkbox"/> pcs 49931	<input type="checkbox"/> pcs 49941	<input type="checkbox"/> pcs 49932
Special variant: mounting bar extended +10cm		<input type="checkbox"/> 12780-1	



16. TRAY TABLE

Info Sheet 95930-1 Tray tables

Tray table, PC (smoke coloured)	<input type="checkbox"/>	sw38 69521	sw42 69531	sw46 69541	sw50 69551	
Tray table, PC (smoke), with pads	<input type="checkbox"/>	69522	69532	69542	69552	
Tray table, Wood (oak)	<input type="checkbox"/>	69525	69535	69545	69555	
Tray table, Wood (oak), with pads	<input type="checkbox"/>	69526	69536	69546	69556	
Support table, padded	<input type="checkbox"/>	----- 69596 -----		----- 69598 -----		
Mug tray (oak)	left <input type="checkbox"/>	69581	right <input type="checkbox"/>	69582		
Tray table mounts, pair	<input type="checkbox"/>	69515-1				
Tray table mounts w/ locking, pair	<input type="checkbox"/>	69515-5				
Mug tray mount w/ locking	left <input type="checkbox"/>	69585-1	right <input type="checkbox"/>	69585-2		
Mounting of tray mounts	<input type="checkbox"/>	69515-8				

17. SPOKE PROTECTION & WHEEL COVERS

Info Sheet 95938-1 Wheels

Spoke protection for model 24, pair	<input type="checkbox"/>	(24" - 15471, 20" - 15472)		
Mounting of spoke protection	<input type="checkbox"/>	15475-8		
Wheel covers for model 16, pair	<input type="checkbox"/>	15690		
Mounting of wheel covers	<input type="checkbox"/>	15690-8		

18. OTHER ACCESSORIES

Profile stop pcs 12799
Assists in mounting supports (head/neck/thoracic) in correct position. Info Sheet 95932.



Bag hook *Info Sheet 95957-1.* pcs 13605



Mounting of bag hook 13605-8

IV pole 69760
Info Sheet 95956-1 and Risk Information 95707-1.



IV pole mount left 69755-1
 right 69755-2



Mounting of IV pole mount 69750-8

Oxygen tank holder 2 litre 59605
Info Sheet 95955-1 and Risk Information 95707-1.



Mounting of oxygen tank holder 56904-8

19. NOTES

