

Why choose HD Balance?

CARER ERGONOMY

Easy to drive and maneuver. Easy to tilt, even with patient that has more substantial body weight. Convenient placement of levers for tilt and back recline.

CHOCK ABSORBING BACK

Possibility to get a chock absorbing back, suitable for example to spastic patients. Spares the patient's body and the wheelchair's frame.

FLEXI BACK

Unique backrest that can adjust the pelvis position and support the lower back. Can also be shaped after Kyphosis, wide bottoms and u-shaped around the back for improved lateral support.

Correct adjustment of the backrest can help improve breathing, swallowing and speech.

ADAPTABLE BACK WIDTH

Possible to get a back width that is one or two sizes wider than the seat width. Also one size narrower is available.

SCALES/USER CARD

Has scales for back recline, seat tilt, armrest height, legrest angle, and seat depth.

User card eases the communication between carers and facilitates a good patient security.



LONG LIFE TIME

Expected to be 10 years at normal use.

DEPTH ADJUSTABLE ARMRESTS

Possibility to put back the armrest (without the use of tools) gives support for the elbows even in reclined position.

DURABLE


Durable design as standard and also available in strengthened edition to meet demands from, for example, spastic patients.

SIMPLE TO WORK WITH

Easy to adjust for the prescriber, most adjustments are without tools or with the tool included under the seat.

Easy to adjust for the carer. The headrest, for example, is adjusted for the current posture without tools.

Easy to repair and recondition for technicians.



LOW TOTAL WEIGHT

Weighs only 30-31 kg including carer brakes, armrests, angle adjustable legrests and cushions.

ANGLE ADJUSTABLE LEGRESTS 90-180°

Adjustable between 90-180° without interfering with the castors. Also available with fixed legrests, 90° and 105°.

KEPT POSTURE WHILE REMOVING/ ATTACHING LEGRESTS

Legrests can be removed/attached behind the patients calves without straightening the legs and change the posture.

Particularly good for patients with stretch spasticity or contractures.

ADAPTABLE

HD Balance is available in many editions and sizes and a very extensive range of accessories.

Check our web page for more info:

www.hdrehab.se